

ARCADIA MEDICAL®

Silicone Cuffless Extended Connect® Pediatric & Neonatal Tracheostomy Tubes




Latex Free

Exceptional Comfort In Airway Management

ARCADIA MEDICAL®

CE
0297

Silicone Cuffless Extended Connect® Pediatric & Neonatal Tracheostomy Tubes



Arcadia Medical® Silicone Cuffless Extended Connect® Pediatric and Neonatal Tracheostomy Tubes are wire reinforced along the entire tube length providing exceptional comfort and safe, kink resistant airway management.

Extended Connect® extends the connection away from the chin to help prevent abrasion, enhancing patient comfort and mobility. Extended Connect® also protects against tube occlusion associated with flexing of the head and neck.

Soft, flexible, biocompatible silicone is the ideal material to help protect the sensitive airway of Pediatric and Neonatal patients.

Silicone is hydrophobic, naturally lubricious and less susceptible to encrustation. Silicone is latex and plasticizer free.

Arcadia Medical® Silicone Cuffless Extended Connect® Pediatric Tracheostomy Tubes (Wire Reinforced)

Reference Number	Tube Size	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Tube Length (mm) (Machine End)	Distal Tube Length (mm) (Patient End)	Units Per Box
152 025	2.5	2.5	4.0	20	38	1
152 030	3.0	3.0	4.7	20	39	1
152 035	3.5	3.5	5.3	25	40	1
152 040	4.0	4.0	6.0	25	41	1
152 045	4.5	4.5	6.7	30	42	1
152 050	5.0	5.0	7.3	30	44	1
152 055	5.5	5.5	8.0	30	46	1

Arcadia Medical® Silicone Cuffless Extended Connect® Neonatal Tracheostomy Tubes (Wire Reinforced)

Reference Number	Tube Size	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Tube Length (mm) (Machine End)	Distal Tube Length (mm) (Patient End)	Units Per Box
157 025	2.5	2.5	4.0	15	30	1
157 030	3.0	3.0	4.7	17	32	1
157 035	3.5	3.5	5.3	20	34	1
157 040	4.0	4.0	6.0	22	36	1

Packaged individually sterile with obturator, twill tie and disconnect wedge.

Made in USA

Latex Free

US AGENT

EC REP

Arcadia Medical Corporation
 1450 E. American Lane, Suite 1400
 Schaumburg, IL 60173 USA
 Phone 847.330.4447
 Fax 847.620.2502

RMS-UK, Ltd.
 28 Trinity Road
 Nailsea, Somerset
 BS48 4NU, UK

E-Mail info@arcadiamedical.com
 Website www.arcadiamedical.com

0297

Corporate Office

Arcadia Medical Corporation
 Craigmuir Chambers, Russell Hill Road
 P.O. Box 71, Road Town, Tortola
 British Virgin Islands