

ARCADIA MEDICAL®

Silicone Cuffless Pediatric & Neonatal Tracheostomy Tubes



Latex Free


Exceptional Comfort In Airway Management

ARCADIA MEDICAL®

CE
0297

ARCADIA MEDICAL®

Silicone Cuffless Pediatric & Neonatal Tracheostomy Tubes



Arcadia Medical® Silicone Cuffless Pediatric and Neonatal Tracheostomy Tubes are wire reinforced to provide exceptional comfort and safe, kink resistant airway management.

Soft, flexible, biocompatible silicone is the ideal material to help protect the sensitive airway of Pediatric and Neonatal patients.

Silicone is hydrophobic, naturally lubricious and less susceptible to encrustation. Silicone is latex and plasticizer free.

The angled soft silicone neck flange is comfortable to wear and minimizes stomal irritation. The neck flange has an integral 15mm ISO Swivel Connector.

Arcadia Medical® Silicone Cuffless Pediatric Tracheostomy Tubes (Wire Reinforced)

Reference Number	Tube Size	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Units Per Box
150 025	2.5	2.5	4.0	38	1
150 030	3.0	3.0	4.7	39	1
150 035	3.5	3.5	5.3	40	1
150 040	4.0	4.0	6.0	41	1
150 045	4.5	4.5	6.7	42	1
150 050	5.0	5.0	7.3	44	1
150 055	5.5	5.5	8.0	46	1

Arcadia Medical® Silicone Cuffless Neonatal Tracheostomy Tubes (Wire Reinforced)

Reference Number	Tube Size	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Units Per Box
155 025	2.5	2.5	4.0	30	1
155 030	3.0	3.0	4.7	32	1
155 035	3.5	3.5	5.3	34	1
155 040	4.0	4.0	6.0	36	1

Packaged individually sterile with obturator, twill tie and disconnect wedge.

Made in USA

US AGENT

Arcadia Medical Corporation

1450 E. American Lane, Suite 1400
Schaumburg, IL 60173 USA
Phone 847.330.4447
Fax 847.620.2502

E-Mail info@arcadiamedical.com
Website www.arcadiamedical.com

Corporate Office

Arcadia Medical Corporation

Craigmuir Chambers, Russell Hill Road
P.O. Box 71, Road Town, Tortola
British Virgin Islands

 Latex Free

EC REP

RMS-UK, Ltd.

28 Trinity Road
Nailsea, Somerset
BS48 4NU, UK


0297

ARCADIA MEDICAL®

1011 Rev. B
Arcadia Medical® is a trademark of
Arcadia Medical Corporation